



Vista Plus

Transilluminator - Easy Identification Of Veins



The Vista Plus is a portable vein illumination device for both paediatric and adult patients, and patients with a history of difficult venous access.

Vista Plus also improves the visualisation of veins in geriatric patients and patients with range of pigmented skin types.

When activated it reveals a map of the vasculature on the surface of the skin. This leads to a more effective venepuncture procedure, with less patient discomfort, resulting in higher patient satisfaction and ensuring shorter procedure time.

- Easy to use, economical and hands free.
- Increases the success rate of first attempt venepuncture.
- 3 illumination settings for maximum patient versatility.
- A rechargeable lithium battery provides 40 min of continuous run time at full illumination and recharges in 2hrs.
- Special indicators show the remaining charge.
- Splash proof, can be wiped over with disinfectant.

Did You Know?

- More than 50% of children admitted to a healthcare facility requires the insertion of a peripheral intravenous cannula (PIVC).
- More than 50% of these insertions fail on the first attempt.

SIZES AVAILABLE

When ordering please quote the following (minimum order qty 1 box):

ORDER#	DIMENSIONS	QTY/box
TRANS2	162mm x 40mm x 35mm	1

Distributed By:



T: 1300 106 003

07 3117 9539

F: 07 3036 6625

E: sales@hmsmedical.com.au